# DiaMetra

## Anti-Tg

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
## Anti-Tg


for routine analysis

Quantitative determination of antibodies against Thyroglobulin in human serum or plasma.

LOT

See external label

Σ = 96 tests

REF DKO115

### INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Anti-Tg concentration in serum and plasma.

### 1. CLINICAL SIGNIFICANCE

Thyroglobulin (Tg) is a protein produced by and used entirely within the thyroid gland. Thyroglobulin is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Thyroglobulin synthesis is stimulated at the transcriptional level by thyroid-stimulating hormone (TSH). Thyroglobulin (Tg) is a well-known target for autoantibodies occurring in thyroid autoimmunity (Graves' disease and Hashimoto's thyroiditis). Anti-Tg antibodies mostly belong to the IgG class. Low to moderate levels of anti-Tg antibodies can be found in sera of other autoimmune patients (eg systemic lupus erythematosus or Sjogren syndrome). In some cases anti-Tg positive sera may show negativity for other type of anti-thyroid antibodies - anti-TPO.

Therefore, combined determination of both types of anti-thyroid antibodies (anti-TPO + anti-Tg) provides most sensitive laboratory diagnostic tool for thyroid autoimmunity. Separately from autoimmunity, anti-Tg antibodies may develop in patients suffering from thyroid cancer. High level of anti-Tg in such patients may interfere with correct determination of serum thyroglobulin, which serves as tumour marker for therapy control in this group of patients.

### 2. PRINCIPLE

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards human Ig and labeled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of specific antibodies in the specimen.

### 3. REAGENT, MATERIAL AND INSTRUMENTATION

#### 3.1. Reagent and material supplied in the kit

9. Anti-Tg Standards 6x (1 vial = 1.2 mL)
- |      |                    |
|------|--------------------|
| STD0 | REF DCE002/11506-0 |
| STD1 | REF DCE002/11507-0 |
| STD2 | REF DCE002/11508-0 |
| STD3 | REF DCE002/11509-0 |
| STD4 | REF DCE002/11510-0 |
| STD5 | REF DCE002/11511-0 |
10. Controls 2x (1 vial = 1.2 mL)
- |                  |                    |
|------------------|--------------------|
| Negative Control | REF DCE045/11501-0 |
| Positive Control | REF DCE045/11502-0 |
11. Sample Diluent (1 bottle) 100 mL  
Phosphate buffer 50 mM pH 7.4; BSA 1 gr/L  
REF DCE053-0
12. Conjugate (1 bottle) 15 mL  
Human IgG-HRP  
REF DCE002/11502-0
13. Coated Microplate  
(1 microplate breakable)  
REF DCE002/11503-0
14. TMB-Substrate (1 bottle) 15 mL  
H<sub>2</sub>O<sub>2</sub>-TMB (0,25 gr/L) (avoid any skin contact)  
REF DCE004/-0
15. Concentrated Wash Solution (1 bottle) 20 mL  
NaCl 9 gr/l; Tween20 1gr/l  
REF DCE006-0
16. Stop solution (1 bottle) 15 mL (avoid any skin contact)  
0.15 M Sulfuric acid  
REF DCE005-0

#### 3.2. Reagents necessary not supplied

Distilled water.

#### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.  
Microplates reader

### 4. PRECAUTION

- This kit is intended for in vitro use by professional persons only.
- Please adhere strictly to the sequence of pipetting steps provided in this protocol.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- All reagents should be stored refrigerated at 2 - 8 °C in their original container.
- Do not interchange kit components from different lots. The expiration dates stated on the labels of the shipping container and all vials have to be observed. Do not use kit components beyond their expiration dates.
- Allow all kit components and specimen to reach room temperature prior to use and mix well.
- All human source material used in the preparation of standards and controls for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV by FDA cleared methods. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Positive and Negative Control should be handled in the same manner as potentially infectious material.
- Avoid contact between the buffered Peroxide Solution and easily oxidized materials; extreme temperatures may initiate spontaneous combustion.

## 5. PROCEDURE

### 5.1. Preparation of the Standard (S<sub>0</sub>,...,S<sub>5</sub>)

The assay system is calibrated in relative arbitrary units. The standard have approximatively the following concentration:

	S <sub>0</sub>	S <sub>1</sub>	S <sub>2</sub>	S <sub>3</sub>	S <sub>4</sub>	S <sub>5</sub>
AU/mL	0	2	4	8	32	128

### 5.2. Preparation of the Sample

For determination of Anti-Tg serum or plasma are the preferred sample matrixes.

All serum and plasma samples have to be prediluted with sample diluent 1 : 100. Therefore 10 µL of sample may be diluted with 1,000 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum or plasma from the cells by centrifugation after clot formation.

Samples may be stored refrigerated at 2 - 8 °C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20 °C. To avoid repeated thawing and freezing the samples should be aliquoted.

Do not use microbiologically contaminated samples, as highly lipemic or hemolized samples.

### 5.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50x) with distilled water to a final volume of 1000 ml prior to use. Store refrigerated: stable at 2 - 8 °C for at least 30 days after preparation or until the expiration date printed on the label.

### 5.4. PROCEDURE

Allow all reagents to stand at room temperature (22-28°C). As it is necessary to perform the determination in duplicate, prepare two wells for each of the six points of the standard curve (S<sub>0</sub>-S<sub>5</sub>), two for each sample and control, one for Blank.

Reagent	Standard	Sample	Blank
Standard S <sub>0</sub> -S <sub>5</sub>	100 µL		
Controls	100 µL		
Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL diluted wash solution			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL diluted wash solution.			
TMB substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop solution	100 µL	100 µL	100 µL
Read the absorbance (E) at 450 nm against Blank.			

## 6. RESULTS

### 6.1. Standard curve

For Anti-Tg a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable.

### Recommended Lin-Log Plot

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

### TYPICAL RESULTS (example only)

The figures below show typical results for Anti- Tg. These data are intended for illustration only and should not be used to calculate results from another run.

STD	OD1	OD2	mean	C1	C2	mean	CV%
0	0.025	0.032	0.0285	0.045	0	0.0225	141,42
1	0.231	0.24	0.2355	1.88	2.12	2	8,49
2	0.429	0.442	0.4355	3.96	4.2	4.08	4,16
3	0.767	0.759	0.763	8.2	7.77	7.985	3,81
4	1.884	2.065	1.9745	31.9	32.36	32.13	1,01
5	3.190	3.316	3.253	128.2	127.5	127.85	0,39

**DiaMetra S.r.l. Headquater:** Via Garibaldi, 18 – 20090 SEGRATE (MI)

Tel. 0039-02-2139184 – 02-26921595

Fax 0039-02-2133354.

**Manufact:** Via Giustozzi, 35/35a (già via Bartolomei) – Z.I Paciana – 06034 FOLIGNO (PG) ITALY.

Tel. 0039-0742-24851

Fax 0039-0742-316197

E-mail: info@diametra.com

## 7. REFERENCE VALUE

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-Tg tests:

Anti-Tg [AU/ml]

normal: < 4

elevated: ≥ 4

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-Tg.

### 7.1. Specificity

Comparison test against a commercial reference kit, assayed on 66 serums (26 of them positive serums and 40 negative serums) shows a 95% specificity.

### 7.2. Sensibility

Comparison test against a commercial reference kit, assayed on 66 serums (26 of them positive serums and 40 negative serums) shows a 88.5% sensibility.

### 7.3. Detection limit

The lowest concentration of anti-Tg that can be distinguished from zero standard is about 0.31 AU/mL with a confidence limit of 98%.

## 8. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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### **ERROR POSSIBLE CAUSES / SUGGESTIONS**

#### **No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

#### **Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

#### **Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

#### **Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation