DiaMetra

Anti-TPO

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
Anti-TPO

for routine analysis

Quantitative determination of antibodies against Thyroperoxidase in human serum or plasma.

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO116

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Anti-TPO concentration in serum and plasma.

1. CLINICAL SIGNIFICANCE

Anti-TPO antibodies (formerly thyroid microsomal antibodies) are directed against a target protein - thyroid peroxidase (TPO) - located in the smooth endoplasmic reticulum of thyroid cells.

Thyroid peroxidase or Thyroperoxidase (TPO) is an enzyme mainly expressed in the thyroid that liberates iodine for addition onto tyrosine residues on thyroglobulin for the production of thyroxine (T4) or triiodothyronine (T3) (thyroid hormones).

It is a frequent epitope of autoantibodies in autoimmune thyroid disease. In autoimmune disorders, your immune system makes antibodies that mistakenly attack normal tissue such as thyroid peroxidase. Antibodies that attack the thyroid gland cause inflammation and impaired function of the thyroid.

The presence of anti-TPO antibodies in serum is associated with thyroid autoimmune diseases (Graves' disease and Hashimoto's thyroiditis). Anti-TPO antibodies mostly belong to the IgG class.

Low to moderate levels of serum anti-TPO antibodies can be found in some other autoimmune pathology (eg systemic lupus erythematosus or Sjogren syndrome) and, rarely, in apparently healthy subjects (especially elderly women). Anti-TPO antibodies are more sensitive in diagnostics of thyroid autoimmune diseases than anti-thyroglobulin (anti-TG) antibodies.

However, in some cases anti-TG positive sera may be negative for anti-TPO. Therefore, combined determination of both types of anti-thyroid antibodies (anti-TPO + anti-TG) provides a more sensitive laboratory diagnostic tool for thyroid autoimmunity.

2. PRINCIPLE

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards human Ig and labeled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen substrate mixture, stop solution and photometry

at 450 nm. Optical density in the microwell is directly related to the quantity of specific antibodies in the specimen.

3. REAGENT, MATERIAL AND INSTRUMENTATION

3.1. Reagent and material supplied in the kit

1. Anti-TPO Standards 6x (1 vial = 1.2 mL)
Phosphate buffer 0.1M, NaN₃ <0.1%

STD0	REF DCE002/11606-0
STD1	REF DCE002/11607-0
STD2	REF DCE002/11608-0
STD3	REF DCE002/11609-0
STD4	REF DCE002/11610-0
STD5	REF DCE002/11611-0

2. Controls 2x (1 vial = 1.2 mL)
Phosphate buffer 0.1M, NaN₃ <0.1%

<u>Negative Control</u>	REF DCE045/11601-0
<u>Positive Control</u>	REF DCE045/11602-0

3. Sample Diluent (1 bottle) 100 mL
Phosphate buffer 0.1M, NaN₃ <0.1%
REF DCE053-0

4. Conjugate (1 bottle) 15 mL
Anti h-IgG conjugated with peroxidase, BSA 0.1%, Proclin <0.003%
REF DCE002/11602-0

5. Coated Microplate
(1 breakable microplate) Thyroperoxidase adsorbed on microwells surface
REF DCE002/11603-0

6. TMB-Substrate (1 bottle) 15 mL
3,3',5,5'-tetramethylbenzidine 0.26 g/L,
hydrogen peroxide 0.05%
REF DCE004-0

7. Concentrated Wash Solution (1 bottle) 20 mL
Phosphate buffer 0.2M, proclin <0.002%
REF DCE006-0

8. Stop solution (1 bottle) 15 mL
Sulfuric acid 0.15M
REF DCE005-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.
Microplates reader

4. PRECAUTION

- This kit is intended for in vitro use by professional persons only.
- Please adhere strictly to the sequence of pipetting steps provided in this protocol.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- All reagents should be stored refrigerated at 2 - 8 °C in their original container.
- Do not interchange kit components from different lots. The expiration dates stated on the labels of the shipping container and all vials have to be observed. Do not use kit components beyond their expiration dates.
- Allow all kit components and specimen to reach room temperature prior to use and mix well.
- All human source material used in the preparation of standards and controls for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV by FDA cleared methods. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Positive and Negative Control should be handled in the same manner as potentially infectious material.
- Avoid contact between the buffered Peroxide Solution and easily oxidized materials; extreme temperatures may initiate spontaneous combustion.

5. STORAGE CONDITION

- Store all the kit reagents at 2-8°C. Do not freeze. Reagents are stable until the expiration date when stored and handled as directed.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.

6. PROCEDURE

6.1. Preparation of the Standard (S₀,...,S₅)

The assay system is calibrated in relative arbitrary units. The standard have approximatively the following concentration:

	S ₀	S ₁	S ₂	S ₃	S ₄	S ₅
AU/mL	0	5	10	20	80	320

6.2. Preparation of the Sample

For determination of Anti-TPO serum or plasma are the preferred sample matrixes.

All serum and plasma samples have to be prediluted with sample diluent 1 : 100. Therefore 10 µL of sample may be diluted with 1,000 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum or plasma from the cells by centrifugation after clot formation.

Samples may be stored refrigerated at 2 - 8 °C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20 °C. To avoid repeated thawing and freezing the samples should be aliquoted. Do not use microbiologically contaminated samples, as highly lypemic or hemolized samples.

6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50x) with distilled water to a final volume of 1000 ml prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8 °C.

6.4. PROCEDURE

Allow all reagents to stand at room temperature (22-28°C). As it is necessary to perform the determination in duplicate, prepare two wells for each of the six points of the standard curve (S₀-S₅), two for each sample and control, one for Blank.

Reagent	Standard	Sample	Blank
Standard S ₀ -S ₅	100 µL		
Controls	100 µL		
Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL diluted wash solution			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL diluted wash solution.			
TMB substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop solution	100 µL	100 µL	100 µL
Read the absorbance (E) at 450 nm against Blank.			

7. RESULTS

7.1. Standard curve

For Anti-TPO a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction

method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable.

Recommended Lin-Log Plot

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

TYPICAL RESULTS (example only)

The figures below show typical results for Anti-TPO. These data are intended for illustration only and should not be used to calculate results from another run.

STD	OD1	OD2	mean	C1	C2	Mean	CV%
0	0,012	0,013	0,013	0	0,03	0,02	141,42
1	0,125	0,121	0,123	5,65	5,45	5,55	2,49
2	0,212	0,210	0,211	9,94	9,84	9,89	0,71
3	0,395	0,422	0,409	19,37	20,83	20,10	5,12
4	1,303	1,221	1,262	84,12	76,43	80,28	6,77
5	2,521	2,474	2,498	330,2	310,9	320,52	4,26

8. REFERENCE VALUE

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-TPO tests:

Anti-TPO [AU/ml]

normal: < 20

elevated: ≥ 20

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-TPO.

8.1. Specificity

Comparison test against a commercial reference kit, assayed on 78 serums (48 of them positive serums and 30 negative serums) shows a 96.6% specificity.

8.2. Sensibility

Comparison test against a commercial reference kit, assayed on 78 serums (48 of them positive serums and 30 negative serums) shows a 97.9% sensibility.

8.3. Detection limit

The lowest concentration of anti-TPO that can be distinguished from zero standard is about 0.31 AU/mL with a confidence limit of 95%.

9. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation