# Free PSA

For use on IMMULITE® 2000 systems

**SIEMENS**

# IMMULITE® 2000 Free PSA

## English

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of free, uncomplexed prostate-specific antigen (PSA), i.e. PSA not bound to  $\alpha_1$ -antichymotrypsin or other binding proteins, in serum.

Catalog Number: **L2KPF2** (200 tests)

Test Code: **fPS** Color: **Light Gray**

The concentration of free PSA in a given specimen determined with different assays can vary due to the differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the assay used. Values obtained with different PSA assays can not be used interchangeably. Before changing assays, the laboratory must confirm the baseline values for patients being serially monitored.

## Summary and Explanation

Since the inception of prostate-specific antigen (PSA) measurements by immunoassay, it has been observed that several distinct immunoreactive species could be isolated from the serum of prostate cancer (PCa) patients after molecular sieve chromatography.<sup>1</sup> PSA has been shown to form stable complexes with two of the major extracellular protease inhibitors in blood,  $\alpha_1$ -antichymotrypsin (ACT) and  $\alpha_2$ -macroglobulin (AMG).<sup>4</sup> In prostate cancer patients, PSA complexed with ACT (PSA-ACT) is typically the major form in circulation: for about 50% of these patients, PSA-ACT accounts for 85% of the total PSA present. Some 12–15% of prostate cancer patients, on the other hand, present with free (that is, uncomplexed) PSA as the predominant form.<sup>2</sup>

Biochemical studies of PSA isolated from seminal plasma show that approximately

35% is free, being enzymatically inactive and unreactive with protease inhibitors.<sup>3</sup> According to current hypotheses, this form of the molecule represents *either* an inactive zymogen of PSA *or* a nicked or enzymatically inactive form.

Complex formation with ACT results in exposure of a limited number of the antigenic epitopes of PSA, whereas complex formation with AMG encapsulates the antigenic epitopes of PSA. Differences in recognition of these multiple forms of PSA by reagent antibodies have contributed to discrepancies between commercial PSA assays.<sup>4</sup>

Immunometric assays have now been designed to selectively characterize all the molecular forms of PSA in circulation: some detect just PSA-ACT, others just free PSA, and some are considered to be assays for total PSA, detecting PSA epitopes available both on free (uncomplexed) PSA and on PSA complexed to serine protease inhibitors.

Using assays such as these, free PSA was found to comprise a significantly ( $p < 0.0001$ ) smaller fraction of the total PSA concentration in patients with untreated PCa than in patients with benign prostatic hypertrophy (BPH).<sup>2</sup>

A substantial body of literature indicates the benefit of combining free and total PSA in the form of a ratio, facilitating discrimination between PCa and BPH. These determinations must be made using assays developed by the same manufacturer.<sup>14,15</sup> The free-to-total ratio (% free PSA) is typically expressed as a percent:  $100 \times \text{free PSA} / \text{total PSA}$ , with both measured in ng/mL. The % free PSA is, on average, lower in PCa than in BPH and has been commonly used as an aid in PCa diagnosis when the total PSA concentration falls in the “gray zone”, that is, between 4 and 10 ng/mL.<sup>16–19</sup> While a lower cutoff might result in fewer false-positive results, a higher cutoff is less likely to miss actual instances of PCa. The % free PSA has also been investigated in the wider total PSA concentration range of 2.0 to 20 ng/mL.<sup>15,17,18,20</sup>

The optimal % free PSA cutoff depends on other factors as well, one of which is the prostate and/or transition zone volume measured by transrectal ultrasound (TRUS).<sup>21,22</sup> Age is another factor: while both free and total PSA rise with advancing age, the % free PSA decreases.<sup>23</sup> More complex approaches employing artificial neural networks (ANNs) take these parameters and others into account, e.g., the results of a digital rectal exam (DRE). Using an ANN can significantly increase discrimination between BPH and PCa.<sup>24</sup>

## Principle of the Procedure

IMMULITE 2000 Free PSA is a solid-phase, sequential chemiluminescent immunometric assay.

The solid phase (bead) is coated with monoclonal murine anti-PSA antibody specific for free (uncomplexed) PSA. The patient sample is incubated with the bead during the first cycle at which time free PSA in the sample binds to the free PSA-specific monoclonal antibody-coated bead. Unbound serum is then removed by a centrifugal wash. Alkaline phosphatase (bovine calf intestine) conjugated to a goat anti-PSA polyclonal antibody in the reagent is introduced in the second cycle and binds to free PSA on the bead to form an antibody-sandwich complex. Unbound enzyme conjugate is then removed by a centrifugal wash. Finally, chemiluminescent substrate is added to the bead and signal is generated in proportion to the bound enzyme.

**Incubation Cycles:** 2 × 30 minutes

**Time to first result:** 65 minutes

## Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Free PSA has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 25 µL serum

**Storage:** 24 hours at 2–8°C or for longer at –20°C.<sup>29,30</sup>

## Warnings and Precautions

For *in vitro* diagnostic use only.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>26-28</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H302 + H312,  
H412**

**P280, P273,  
P301 + P312,  
P302 + P312,  
P501**

**Warning!** Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.  
**Contains:** sodium azide; Free PSA Adjustors

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

### Free PSA Bead Pack (L2PF12)

With barcode. 200 beads, coated with monoclonal murine anti-PSA antibody specific for free (uncomplexed) PSA. Stable at 2–8°C until expiration date.

**L2KPF2:** 1 pack

### Free PSA Reagent Wedge (L2PFA2)

With barcodes. 11.5 mL of a protein buffer matrix, with preservative; 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-PSA antibody specific for PSA in a buffer containing human serum, with preservative. Stable at 2–8°C until expiration date.

**L2KPF2:** 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Free PSA Adjustors (LPFL, LPFH)

Two vials (Low and High) of lyophilized free PSA in a buffer solution, with preservative. Reconstitute each vial with **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KPF2:** 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

**L2M2Z:** 25 mL    **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M2Z:** 3 labels    **L2M2Z4:** 5 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also Required  
Distilled or deionized water; test tubes;  
controls

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of free PSA.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

Based on its relationship to IMMULITE Free PSA (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

A reference range study for IMMULITE Free PSA was performed using serum samples from adult volunteers, including men ranging from approximately 20 to 70 years of age (central 95%: 22–64 years, median: 40 years). The subjects were in apparent good health, based on a questionnaire.

Blood samples were collected in France, Germany, the Netherlands and Portugal. Results were generated by an independent laboratory in the Netherlands using IMMULITE kits. (The samples were collected in plain glass tubes without anticoagulants, gel barriers, or clot-promoting additives, and assayed in singlicate.)

Medians and 95th percentiles are tabulated below. Percentiles were determined nonparametrically.

Group	Median	95%ile	Units	<i>n</i>
Males	0.14	0.42	ng/mL	425

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

Serum PSA concentrations should not be interpreted as absolute evidence for the presence or absence of malignant disease.<sup>9</sup>

Prediction of malignant prostatic disease recurrence should be based on a complete clinical evaluation of the patient, which may also include serial serum PSA determinations.

Samples should be obtained before biopsy, prostatectomy or prostatic massage.<sup>10</sup>

PSA expression may be altered due to hormonal therapy for prostate cancer. Consequently, a low PSA result following a prostatic cancer treatment which includes hormonal therapy may not adequately reflect the presence of residual or recurrent disease.<sup>13</sup>

The accuracy of free PSA assays cannot be evaluated by determining the recovery of free PSA spiked into serum samples, because the formation of complexes between free PSA and the PSA-binding proteins present in normal sera ( $\alpha_1$ -antichymotrypsin and  $\alpha_2$ -macroglobulin) will generally cause erroneously low recoveries.

Some individuals have antibodies to mouse protein which can cause interference in immunoassays that employ antibodies derived from mice. Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy, in particular, may contain human anti-mouse antibodies (HAMA). These specimens may show erroneous results in such assays.<sup>6-8</sup> Therefore, IMMULITE 2000 Free PSA results should be used only in conjunction with results from some other diagnostic procedure and information available from the clinical evaluation of the patient.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**Reportable Range:** 0.07 to 25 ng/mL (WHO NIBSC 1<sup>st</sup> IS 96/668)

**Analytical Sensitivity:** Limit of Blank (highest value expected for a sample with no analyte; calculated as the value lying two standard deviations above that of the lowest calibration): 0.02 ng/mL

Limit of Detection (lowest detectable concentration; determined in accordance with CLSI EP17-A<sup>25</sup>): 0.07 ng/mL

**Functional Sensitivity:** (concentration with 20% coefficient of variation (CV) determined in accordance with CLSI EP17-A<sup>25</sup>): 0.07 ng/mL

**High-dose Hook Effect:**  
None up to 15,118 ng/mL

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Specificity:** IMMULITE 2000 Free PSA is highly specific for prostate-specific antigen, with a particularly low crossreactivity to other naturally occurring compounds that might be present in patient samples. (See "Specificity" table.)

**Bilirubin:** Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Biotin:** Specimens that contain biotin at a concentration of 3500 ng/mL demonstrate a less than or equal to 10% change in results.

**Hemolysis:** Presence of packed red blood cells in concentrations up to 30  $\mu$ L/mL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 15 volunteers into plain, heparinized, EDTA and Becton Dickinson SST<sup>®</sup> vacutainer tubes. All samples were assayed by the IMMULITE 2000 Free PSA procedure. By linear regression:

(Heparin) = 1.04 (Serum) – 0.19 ng/mL  
r = 0.997

(EDTA) = 1.03 (Serum) – 0.19 ng/mL  
r = 0.998

(SST) = 1.08 (Plain Tubes) – 0.28 ng/mL  
r = 0.99

Means:

37 ng/mL (Heparin)  
37 ng/mL (EDTA)  
38 ng/mL (SST)  
36 ng/mL (Serum)

**Method Comparison:** The assay was compared to IMMULITE Free PSA on 184 samples. (Concentration range: approximately 0.07 to 25 ng/mL. See graph.) By linear regression:

(IML 2000) = 1.04 (IML) – 0.08 ng/mL  
r = 0.998

Means:

2.37 ng/mL (IMMULITE 2000)  
2.35 ng/mL (IMMULITE)

## References

- 1) Chu TM, et al. Circulating antibody to prostate antigen in patients with prostatic cancer. *Ann NY Acad Sci* 1983;417: 383-9.
- 2) Lilja H, et al. Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha_1$ -antichymotrypsin. *Clin Chem* 1991;37: 1618-24.
- 3) Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine protease inhibitors. *Eur J Biochem* 1990;194:755-63.
- 4) Zhou AM, et al. Multiple forms of prostate-specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. *Clin Chem* 1993;39:2483-91.
- 5) Jacobs DS, Grady HD, editors. *Laboratory Test Handbook*. 4th ed. Hudson (Cleveland): Lexi-Comp Inc., 1996; 193.
- 6) Primus FJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine antibody for diagnosis and therapy. *Clin Chem* 1988;34:261-4.
- 7) Hansen HJ, et al. Solving the problem of antibody interference in commercial "sandwich"-type immunoassay of carcinoembryonic antigen. *Clin Chem* 1989;35:146-51.
- 8) Schroff RJ, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-85.
- 9) Kuriyama M, et al. Quantitation of prostate-specific antigen in serum by a sensitive enzyme immunoassay. *Cancer Res* 1980;40:4658-62.
- 10) Stamey TA, Yang N, et al. Prostate specific antigen as a serum marker for adenocarcinoma of the prostate. *N Engl J Med* 1987;317:909-16.
- 11) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
- 12) Sensabaugh GF, Crim D. Isolation and characterization of a semen-specific protein from human seminal plasma: a potential new marker for semen identification. *J Forensic Sci* 1978;23:106-15.
- 13) Morgan WR et al. Prostate specific antigen values after radical retropubic prostatectomy for adenocarcinoma of the prostate: impact of adjuvant treatment (hormonal and radiation). *J Urol* 1991;145:319-23.
- 14) Wolff JM, Stocker G, Borchers H, Haubeck H, Greiling H, Jakse G. Critical aspects related to the interpretation of the free-to-total PSA-ratio. *Anticancer Res* 1999;19:2633-6.
- 15) Patel D, White PA, Milford Ward A. A comparison of six commercial assays for total and free prostate specific antigen (PSA): the predictive value of the ratio of free to total PSA. *BJU Int* 2000;85:686-9.
- 16) Correale M, Pagliarulo A, Donatuti G, Sturda F, Capobianco AM, Stigliani V, et al. Preliminary clinical evaluation of free/total PSA ratio by the IMMULITE system. *Int J Biol Markers* 1996;11:24-8.
- 17) Wymenga LF, Duisterwinkel FJ, Groenier K, Visser-van Brummen P, Marrink J, Mensink HJ. Clinical implications of free-to-total immunoreactive prostate-specific antigen ratios. *Scand J Urol Nephrol* 2000;34:181-7.
- 18) McArdle PA, Pollock MA, Wallace AM, McMillan DC, Crooks JE, Underwood MA. Comparison of total, complexed and free prostate-specific antigens and their ratios in the detection of prostate cancer in a non-screened population. *Ann Clin Biochem* 2004;41:201-6.
- 19) Martinez-Pineiro L, Garcia Mediero JM, Gonzalez Gancedo P, Tabernero A, Lozano D, et al. Probability of prostate cancer as a function of the percentage of free prostate-specific antigen in patients with a non-suspicious rectal examination and total prostate-specific antigen of 4-10 ng/ml. *World J Urol* 2004;22:124-31.
- 20) Wesseling S, Stephan C, Semjonow A, Lein M, Brux B, Sinha P, et al. Determination of non-alpha1-antichymotrypsin-complexed prostate-specific antigen as an indirect measurement of free prostate-specific antigen: analytical performance and diagnostic accuracy. *Clin Chem* 2003;49:887-94.
- 21) Stephan C, Lein M, Jung K, Schnorr D, Loening SA. The influence of prostate volume on the ratio of free to total prostate specific antigen in serum of patients with prostate carcinoma and benign prostate hyperplasia. *Cancer* 1997 Jan;79:104-9.
- 22) Stamey TA, Yemoto CE. Examination of the 3 molecular forms of serum prostate specific antigen for distinguishing negative from positive biopsy: relationship to transition zone volume. *J Urol* 2000;163:119-26.
- 23) Wolff JM, Borchers H, Rohde D, Jakse G. Age related changes of free and total prostate specific antigen in serum. *Anticancer Res* 1999;19:2629-32.
- 24) Stephan C, Cammann H, Semjonow A, Diamandis EP, Wymenga LF, Lein M, et al. Multicenter evaluation of an artificial neural network to increase the prostate cancer detection rate and reduce unnecessary biopsies. *Clin Chem* 2002;48:1279-87.
- 25) CLSI. *Protocols for the Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI document EP17-A Vol 24 (No 34). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2004.
- 26) Centers for Disease Control. *Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings*. *MMWR*, 1988;37:377-82, 387-8.
- 27) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 28) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*, 29 CFR 1910.1030.
- 29) Woodrum D, French C, Shamel LB. Stability of free prostate-specific antigen in serum samples under a variety of sample collection and sample storage conditions. *Urology* 1996;48(6A):33-39.
- 30) Woodrum D, York L. Two-year stability of

free and total PSA in frozen serum samples.  
Urology 1998;52(2):247-251.

## Technical Assistance

Available outside the United States only.  
For technical assistance, contact your  
National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of  
Siemens Healthcare Diagnostics Products Ltd.  
is certified to ISO 13485.

## Tables and Graphs

### Precision (ng/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	0.83	0.046	5.5%	0.051	6.1%
2	2.0	0.096	4.8%	0.11	5.5%
3	5.9	0.22	3.7%	0.24	4.1%
4	9.3	0.30	3.2%	0.31	3.3%
5	21	0.87	4.1%	0.91	4.3%

### Specificity


Compound <sup>1</sup>	Amount Added <sup>2</sup>	Percent Cross-reactivity <sup>3</sup>
AFP	10,000 ng/mL	ND
CEA	100 ng/mL	ND
Ferritin	10,000 ng/mL	ND
HCG	100,000 mIU/mL	ND
Lactalbumin	10,000,000 ng/mL	ND
PAP	1000 ng/mL	ND
Prolactin	200 ng/mL	ND

ND: not detectable<sup>4</sup>

### Linearity (ng/mL)

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	4 in 4	0.36	—	—
	2 in 4	0.17	0.18	94%
	1 in 4	0.09	0.09	100%
2	16 in 16 <sup>5</sup>	2.9	—	—
	8 in 16	1.5	1.5	100%
	4 in 16	0.72	0.73	99%
	2 in 16	0.33	0.36	92%
	1 in 16	0.15	0.18	83%
3	16 in 16	5.1	—	—
	8 in 16	2.6	2.6	100%
	4 in 16	1.3	1.3	100%
	2 in 16	0.61	0.64	95%
	1 in 16	0.29	0.32	91%
4	16 in 16	8.3	—	—
	8 in 16	4.1	4.2	98%
	4 in 16	2.1	2.1	100%
	2 in 16	1.1	1.0	110%
5	16 in 16	13	—	—
	8 in 16	6.4	6.5	99%
	4 in 16	3.2	3.3	97%
	2 in 16	1.6	1.6	100%
	1 in 16	0.81	0.81	100%
6	16 in 16	20	—	—
	8 in 16	9.7	10.0	97%
	4 in 16	5.0	5.0	100%
	2 in 16	2.5	2.5	100%
1 in 16	1.3	1.3	100%	

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2020 Siemens Healthcare Diagnostics. All rights reserved.
**CONTROL AB** En: Control Antibody