IMMULITE®

FSH

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 FSH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of follicle stimulating hormone (FSH, follitropin) in serum, as an aid in the diagnosis and treatment of pituitary and gonadal disorders.

Catalog Number: **LKFS1** (100 tests)

Test Code: **FSH** Color: **Light Gray**

Summary and Explanation

Follicle stimulating hormone (follitropin, FSH) is secreted by the β -cells of the anterior pituitary under the control of the gonadotropin releasing hormone produced in the hypothalamus. FSH facilitates the development and maintenance of gonadal tissues, which synthesize and secrete steroid hormones. Circulating levels of FSH are controlled by a negative feedback mechanism on the hypothalamus by steroidal hormones. Although FSH and LH are required for normal sexual function in both males and females, the secretory patterns are very different for the two sexes.

In mature females, FSH initiates the growth and development of ovarian follicles. During ovulation, when the follicle is ruptured, the follicle, now called the corpus luteum, secretes estradiol and progesterone, which control the circulating levels of FSH by a negative feedback effect on the hypothalamus. In menopause, with diminished ovarian function, there is a resulting decrease in estradiol secretion. Due to the lack of a negative feedback effect, with diminished estradiol, the circulating FSH levels become significantly increased.

In the mature male, FSH is associated with the stimulation and maintenance of spermatogenesis. Testosterone and estradiol have the role of providing the negative feedback signal to the hypothalamus for controlling the release of FSH. Infertility in males may be due to hypogonadism as a result of primary testicular failure. Testicular failure may be

functional failure to mature, or a result of germ cell damage. Whatever the etiology, the conditions of hypogonadism have the net result of dramatically raising the circulating FSH levels, due to the lack of a negative feedback effect.

Principle of the Procedure

IMMULITE/IMMULITE 1000 FSH is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma should not be used as a sample type.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 FSH has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 μ L serum (Sample cup must contain at least 100 μ L more than the total volume required.)

Storage: 7 days at 2–8°C or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁵⁻¹⁷

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

FSH Test Units (LFS1)

Each barcode-labeled unit contains one bead coated with monoclonal murine anti-FSH. Stable at 2–8°C until expiration date.

LKFS1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

FSH Reagent Wedge (LFS2)

With barcode. 7.5 mL alkaline phosphatase (bovine calf intestine) conjugated to murine monoclonal anti-FSH antibody in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKFS1: 1 wedge

FSH Adjustors (LFSL, LFSH)

Two vials (Low and High), 3.0 mL each, of FSH in a nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKFS1: 1 set

Kit Components Supplied Separately

FSH Sample Diluent (LFSZ)

For the manual dilution of patient samples. One vial containing 25 mL of an FSH-free nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

Also Required

Sample transfer pipets, distilled or deionized water, controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of FSH.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Reference ranges were generated using IMMULITE FSH (mono/poly) in a multi-national study involving women in apparent good health (age: 16–44 years), who volunteered to have blood samples drawn, on a daily basis, throughout one complete ovulatory cycle. (See "Menstrual Cycle" graph.)

Ovulatory Cycles	n*	FSH, mIU/mL	
		Median	Central 95%
Follicular Phase	54 (762)	6.2	2.8–11.3
Follicular Phase, Days 2 to 3	54 (108)	6.6	3.0–14.4
Midcycle	54 (54)	13.6	5.8–21
Luteal Phase	54 (604)	3.4	1.2–9.0

*Number of subjects (total number of results)

The table below summarizes results for adult males and (untreated) postmenopausal females from studies with IMMULITE FSH (mono/poly). The ranges for women on oral contraceptives and for postmenopausal women on ERT are based on the assay's relationship to Coat-A-Count® FSH IRMA.

Group	n	FSH, mIU/mL	
		Median	Central 95%
Adult Males	135	3.8	0.7–11.1
Adult Females:			
Postmenopausal*	76	90.5	21.7–153
Postmenopausal (ERT)	16	27	9.7–111
Oral Contraceptives	12	1.7	ND–4.9

*Preliminary
ND: not detectable

A cross-sectional study of pediatric values performed with IMMULITE FSH (mono/poly) at a "wellness" clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	FSH, mIU/mL	
			Median	Central 95%
Females	Cord	30	ND	
	0.1–3	57	2.3	0.11–13
	4–9	28	0.8	0.11–1.6
Males	Cord	37	0.24	ND–1.2
	0.1–3	72	0.6	ND–5.5
	4–9	31	0.23	ND–1.9
Combined	Cord	67	0.11	ND–1.1
	0.1–3	129	1.1	ND–10
	4–9	59	0.5	ND–1.8

ND: not detectable

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

In certain cases of infertility, treatment with human gonadotropins poses a potential problem for the accurate measurement of FSH levels. The FSH that is administered can cause the patient to produce antibodies to FSH which will interfere directly with the assay.

Because of pulsatile secretion, samples obtained within the same day from the same patient may fluctuate widely within the reference range, reflecting physiological variation rather than errors in technique or methodology.

Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: up to 170 mIU/mL (8.9 ng/mL) [WHO 2nd IRP 78/549] [interim replacement code 94/632]

Analytical Sensitivity: 0.1 mIU/mL

High-dose Hook Effect: None up to 50,000 mIU/mL

Precision: Samples were assayed in duplicate in 40 runs for a total of 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three FSH solutions (150, 600 and 1800 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for FSH. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 600 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 24 volunteers into plastic Becton Dickinson plain serum, heparinized, EDTA and SST[®] vacutainer tubes. All samples were assayed by the IMMULITE FSH procedure

(Heparin) = 1.10 (Serum) + 0.65 mIU/mL

(SST) = 0.94 (Plain Tubes) + 0.12 mIU/mL

(EDTA) = 0.88 (Serum) + 0.57 mIU/mL

Means:

11.7 mIU/mL (Serum)

13.5 mIU/mL (Heparin)

11.1 mIU/mL (SST)

10.8 mIU/mL (EDTA)

EDTA plasma should not be used as a sample type.

Method Comparison: The IMMULITE FSH mono/mono assay was compared to IMMULITE FSH mono/poly assay on 431 samples. (Concentration range: approximately 0.7 to 130 mIU/mL. See graph.) By linear regression:

(IML m/m) = 0.90 (IML m/p) + 0.35 mIU/mL
r = 0.990

Means:

31.4 mIU/mL (IMMULITE mono/mono)

34.5 mIU/mL (IMMULITE mono/poly)

References

1) Santner S, Santen R, Kulin H, Demers L. A model for validation of radioimmunoassay kit reagents: measurement of follitropin and lutropin in blood and urine. *Clin Chem* 1981;27:1892-5. 2) Odell W, et al. Radioimmunoassay for luteinizing hormone in human plasma or serum. *J Clin Invest* 1967;46:248-55. 3) Davidsohn I, Henry J, editors. *Clinical diagnosis by laboratory methods*. 15th ed. Philadelphia: W.B. Saunders, 1974: 704. 4) Nankin H, et al. Repetitive luteinizing hormone elevations in serum of normal men. *J Clin Endo Metab* 1972;33:558-60. 5) Beitens I, et al. Gonadotropin determinations in timed 3-hour urine collections during the menstrual cycle and LHRH testing. *J Clin Endo Metab* 1976;43:46-55. 6) Chipman J, et al. Interrelationship of plasma and urinary gonadotropins: correlations for 24 hours, for sleep/wake periods, and for 3 hours after luteinizing hormone-releasing hormone stimulation. *Clin Endo Metab* 1981;52:225-30. 7) Kulin H, et al. Integration of pulsatile gonadotropin secretion by timed urinary measurements: an accurate and sensitive 3-hour test. *J Clin Endo Metab* 1975;40:783-9. 8) Kulin H, Santner S. Timed urinary gonadotropin measurements in normal infants, children, and adults, and in patients with disorders of sexual maturation. *J Pediatrics* 1977;90:760-5. 9) Urban M, et al. Comparison of estimates of gonadotropin levels by isolated blood samples, integrated blood concentrations, and timed urinary fractions. *J Clin Endo Metab* 1979;48:732-5. 10) Rebar R, Yen S. In: Dorothy Krieger, editor. *Endocrine rhythms*. New York: Raven Press, 1979. 11) Odell WD, et al. Radioimmunoassay for human follicle-stimulating hormone: physiological studies. *J Clin Invest* 1968;47:2551-62. 12) Sutaria UE, Crooke AC. Selection of patients for treatment with human gonadotrophins. *Int J Fertil* 1971;16:42-6. 13) Babson, AL. The IMMULITE Automated Immunoassay System. *J Clin Immunoassay* 1991;14:83-88. 14) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 15) Centers for Disease Control. *Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings*. *MMWR*, 1988;37:377-82, 387-8. 16) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 17) Federal Occupational Safety and Health Administration, *Bloodborne Pathogens Standard*. 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (mIU/mL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	5.7	0.15	2.6%	0.33	5.8%	
2	19.4	0.59	3.0%	1.04	5.4%	
3	38.5	0.89	2.3%	2.13	5.5%	
4	58.8	2.14	3.6%	3.60	6.1%	
5	87.0	3.18	3.7%	5.81	6.7%	

Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	8.95	—	—
	4 in 8	4.50	4.48	100%
	2 in 8	2.36	2.24	105%
	1 in 8	1.16	1.12	104%
2	8 in 8	13.6	—	—
	4 in 8	7.42	6.80	109%
	2 in 8	3.73	3.40	110%
	1 in 8	1.95	1.70	115%
3	8 in 8	27.8	—	—
	4 in 8	14.9	13.9	107%
	2 in 8	7.67	6.95	110%
	1 in 8	3.81	3.48	109%
4	8 in 8	53.8	—	—
	4 in 8	26.8	26.9	100%
	2 in 8	13.9	13.5	103%
	1 in 8	7.11	6.73	106%
5	8 in 8	68.1	—	—
	4 in 8	37.2	34.1	109%
	2 in 8	21.0	17.0	124%
	1 in 8	11.1	8.51	130%
6	8 in 8	109	—	—
	4 in 8	53.1	54.5	97%
	2 in 8	28.7	27.3	105%
	1 in 8	14.7	13.6	108%

Recovery (mIU/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	1.45	—	—
	A	9.42	8.88	106%
	B	30.8	31.4	98%
	C	98.1	91.4	107%
2	—	4.34	—	—
	A	12.9	11.6	111%
	B	34.6	34.1	101%
	C	107	94.1	114%
3	—	4.88	—	—
	A	12.5	12.1	103%
	B	33.4	34.6	97%
	C	105	94.6	111%
4	—	52.8	—	—
	A	60.0	57.7	104%
	B	83.8	80.2	104%
	C	169	140	121%

Specificity


Compound ¹	Added ²		Apparent ³	% Cross reactivity ⁴
	ng/mL	ng/mL	(ng/mL)	
HCG	2000	215	ND	ND
mIU/mL	200,000	21,538	ND	ND
TSH	10	2.0	ND	ND
uIU/mL	200	40.5	ND	ND
LH	75	7.6	ND	ND
mIU/mL	300	30.5	ND	ND
Prolactin	25	25	ND	ND
ng/mL	200	200	ND	ND
HGH	5.0	5.0	ND	ND
ng/mL	50	50	ND	ND
HPL	10,000	10,000	ND	ND
ng/mL	50,000	50,000	ND	ND

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.
CONTROL AB