IMMULITE®

HCG

For use on the IMMULITE®
and IMMULITE® 1000 systems

SIEMENS

IMMULITE®/IMMULITE 1000 HCG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® and IMMULITE 1000 Analyzers — for the quantitative measurement of human chorionic gonadotropin (HCG) in serum, and for strictly qualitative determinations in urine, as an aid in the detection of pregnancy.

Catalog Number: **LKCG1** (100 tests)

Test Code: **HCG** Color: **Orange**

Summary and Explanation

Human chorionic gonadotropin (HCG) is a two-chain glycoprotein hormone (MW ~37,000) normally found in blood and urine only during pregnancy. It is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy.

According to the literature, circulating HCG typically reaches levels of approximately 2000 mIU/mL one month after conception. A peak level on the order of 100,000 mIU/mL is attained in the third month, after which a gradual decline sets in. Following delivery, the HCG level normally undergoes rapid descent, reaching nonpregnant concentrations (usually less than 5 mIU/mL) some two weeks later.

Ectopic pregnancies and pregnancies terminating in spontaneous abortion tend to have lower than normal circulating HCG levels, while somewhat higher levels are often seen in multiple pregnancies.

Principle of the Procedure

IMMULITE/IMMULITE 1000 HCG is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 42 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufactures may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE/IMMULITE 1000 HCG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 5 µL serum or urine. (Sample cup must contain at least 100 µL more than the total volume required.)

Urine: A freshly voided, first-morning specimen is preferred.²² Centrifuging or filtering (Whatman #1) is recommended to clear cloudy samples.

Dilution of High Serum Samples:

Serum samples expected to have HCG levels above the assay's calibration range (5000 mIU/mL) should be manually diluted with HCG Sample Diluent before assay.

Alternatively, patient samples exceeding 5000 mIU/mL may be diluted on-board 1-in-41 using HCG Sample Diluent Wedges (LCGW) in conjunction with IMMULITE 1000 Windows®.

Storage:

Serum: 7 days at 2–8°C or 2 months at –20°C.¹⁸

Urine: 48 hours at 2–8°C²² or at –20°C until assayed.¹³

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²³⁻²⁵

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: sodium azide; HCG Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

HCG Test Units (LCG1)

Each barcode-labeled unit contains one bead coated with murine monoclonal anti-HCG. Stable at 2–8°C until expiration date.

LKCG1: 100 units

Allow the Test Unit bags to come to room temperature before opening. Open by cutting along the top edge, leaving the ziplock ridge intact. Reseal the bags to protect from moisture.

HCG Reagent Wedge (LCG2)

With barcode. 7.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal ovine anti-HCG, in buffer, with preservative. Store capped and refrigerated: stable at 2–8°C until expiration date. Recommended usage is within 30 days after opening when stored as indicated.

LKCG1: 1 wedge

HCG Adjustors (LCGL, LCGH)

Two vials (Low and High) 2 mL each, of HCG in an HCG-free human serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LKCG1: 1 set

Kit Components Supplied Separately

HCG Sample Diluent (LCGZ, LCGZ4)

For the manual dilution of high patient samples. One vial of ready-to-use, processed, HCG-free human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

LCGZ: 25 mL **LCGZ4:** 100 mL

HCG Sample Diluent Module (L1KCGW1,5)

Required for the on-board dilution of high patient samples using IMMULITE 1000 Windows. Sample Diluent wedges (LCGW) containing HCG-free human serum matrix. Stable at 2–8°C for 90 days after opening.

L1KCGW1: 2 wedges

L1KCGW5: 10 wedges

LSUBX: Chemiluminescent Substrate

LPWS2: Probe Wash Module

LKPM: Probe Cleaning Kit

LCHx-y: Sample Cup Holders (barcoded)

LSCP: Sample Cups (disposable)

LSCC: Sample Cup Caps (optional)

DCHS: Dilution Cup Holders (barcoded)
for IMMULITE 1000 Windows on-board
dilutions

LCGCM: High Level Control Module

Also Required

Sample transfer pipets, distilled or
deionized water, controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE or IMMULITE 1000 Operator's Manual.

See the IMMULITE or IMMULITE 1000 Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Visually inspect each Test Unit for the presence of a bead before loading it onto the system.

Recommended Adjustment Interval:

4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or sample pools with at least two levels (low and high) of HCG.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Qualitative Procedure for Urine

The qualitative procedure allows for determining whether a urine sample has an HCG concentration above or below 30 mIU/mL, within the precision of the assay. Concentrations above this cutoff

may be considered positive for pregnancy providing other causes of high HCG levels have been ruled out, viz. trophoblastic and nontrophoblastic neoplasia.^{4,9,10,15}

Concentrations below this cutoff are to be expected for samples from individuals who are not pregnant, as well as for women in such an early stage of pregnancy that the HCG concentration has not yet reached the decision level set for the procedure. Assaying another sample collected two or more days later should settle the matter, since in early pregnancy HCG levels can be expected to double in approximately this time.¹¹

Interpretation of Results

A result of greater than or equal to 30 mIU/mL indicates that the patient sample is **positive** for pregnancy.

A result of less than 30 mIU/mL indicates that the patient sample is **negative** for pregnancy.

Expected Values

Males and Nonpregnant Females

In a study performed with the IMMULITE HCG assay and involving serum samples from 797 adult volunteers, including 428 males and 369 nonpregnant females, 95% of the results for males were below the assay's detection limit (1.1 mIU/mL) and all were below 2.5 mIU/mL, while 83% of the results for nonpregnant females were below the detection limit, 95% were below 2.7 mIU/mL, and all were below 5.3 mIU/mL. The subjects were in apparent good health based on a questionnaire.

Pregnant Females

A total of 593 serum samples from apparently healthy pregnant women were processed by the IMMULITE HCG assay. The results are summarized below (in mIU/mL) by gestational age (Gest.) – equivalently, by weeks since last menstrual period (LMP). Also see "Expected Values During Pregnancy" graph.

Weeks of Pregnancy

Gest.	LMP	n	Median	Central 95%
1.3–2	3.3–4	30	71	16–156
2–3	4–5	54	607	101–4870
3–4	5–6	34	5243	1110–31,500
4–5	6–7	34	26,983	2560–82,300
5–6	7–8	36	52,090	23,100–151,000
6–7	8–9	33	93,598	27,300–233,000
7–11	9–13	116	117,678	20,900–291,000
11–16	13–18	72	40,989	6140–103,000
16–21	18–23	80	20,868	4720–80,100
21–39	23–41	104	15,352	2700–78,100

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Pathological Conditions

Ectopic pregnancies or abnormal intrauterine pregnancies have circulating HCG levels below normal for that gestational period; after an initial rise these levels plateau or fall.^{3,17,18}

Limitations

The procedure for urine samples is designed solely as a screen for pregnancy, using a 30 mIU/mL cutoff. Hence urinary results should be reported in a strictly qualitative manner, e.g. as "positive" or "negative" for pregnancy — concentration values should *not* be reported for this sample type.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: up to 5000 mIU/mL (WHO 3rd IS 75/537) for quantitative determinations in serum.

Analytical Sensitivity: 1.1 mIU/mL

High-dose Hook Effect:
None up to 600,000 mIU/mL

Intraassay Precision (Within-Run):
Statistics were calculated for samples from the results of *n* replicates in a single run. (See "Intraassay Precision" table.)

Interassay Precision (Run-to-Run):
Statistics were calculated for samples assayed in *n* different runs. (See "Interassay Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three HCG solutions (3280, 6830 and 14,904 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The assay is highly specific for HCG, with low crossreactivity to other glycoprotein hormones present in patient samples. (See "Specificity" table.)

For concentrations as high as 1000 ng/mL, the HCG α -subunit (WHO 75/569) showed less than 1% crossreactivity. HCG β -subunit (WHO 75/551), spiked into HCG Sample Diluent, yielded the tabulated results. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 $\mu\text{L}/\text{mL}$ has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 19 volunteers into plain, heparinized and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of HCG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE HCG procedure.

(Heparin) = 0.92 (Serum) + 2.1 mIU/mL
r = 0.992

(SST) = 0.88 (Plain Tubes) + 7.9 mIU/mL
r = 0.990

Means:

553 mIU/mL (Serum)

513 mIU/mL (Heparin)

499 mIU/mL (SST Tubes)

Method Comparison (Serum): The assay was compared to Coat-A-Count[®] HCG IRMA on 54 serum samples. (Concentration range: approximately 13 to 4450 mIU/mL. See graph.) By linear regression:

(IML) = 1.12 (CAC IRMA) – 7.5 mIU/mL
r = 0.992

Means:

530 mIU/mL (IMMULITE)

485 mIU/mL (Coat-A-Count IRMA)

Method Comparison (Urine): A total of 123 urine samples from pregnant women were processed by the IMMULITE HCG assay and by Double Antibody HCG. The samples, some of which were diluted, had HCG values ranging up to approximately 300 mIU/mL. By linear regression:

IML = 1.2 (DAb) – 15 mIU/mL
r = 0.981

In another study, as a check on the appropriateness of the 30 mIU/mL cutoff used in the qualitative procedure, urinary results for nonpregnant women showed a clean separation from the distribution of repeated measurements made on a 30 mIU/mL sample. (Mean + 6SD for the nonpregnant women: 13 mIU/mL. Mean – 6SD for the 30 mIU/mL sample: 22 mIU/mL.)

References

- 1) Bangham DR, Storring PL. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:390.
- 2) Batzer F. Hormonal evaluation of early pregnancy. *Fertil Steril* 1980;34:1–12.
- 3) Braunstein GD, et al. First-trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Amer J Obstet Gynecol* 1978;131:25–32.
- 4) Clayton L, et al. Endocrine aspects of trophoblastic neoplasia. *J Reprod Med* 1981;26:192–9.
- 5) Corson SL, Horwitz CA, Lau HL, Soderstrom R. Early diagnosis of pregnancy: a symposium. *J Reprod Med* 1981;26:149–78.
- 6) DeCherney AH, Minkin MJ, Spangler S. Contemporary management of ectopic pregnancy. *J Reprod Med* 1981;26:519–23.
- 7) Goldstein DP, Berkowitz RS, Kosasa TS. Measurement of human chorionic gonadotropin. *Major Problems in Obstet Gynecol* 1982;14:64–97.
- 8) Hager H, Van Weemen BK. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:629.
- 9) Kohorn EI. Hydatidiform mole and gestational trophoblastic disease in southern Connecticut. *Obstet and Gynecol* 1982;59:78–84. *Erratum:* The first diagram belongs with the caption to Figure 3, and vice versa.
- 10) Kosasa TS. Measurement of human chorionic gonadotropin. *J Reprod Med* 1981;26:201–6.
- 11) Lenton E, Neal L, Sulaiman R. Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773–8.
- 12) Lindstedt G, Janson PO, Thorburn J. Sensitivity of serum chorionic gonadotropin assay for ectopic pregnancy. *Lancet* 1981;1:781–2.
- 13) McCready J, Braunstein GD, et al. Modification of the choriogonadotropin beta-subunit radioimmunoassay for determination of urinary choriogonadotropin. *Clin Chem* 1978;24:1958–61.
- 14) Storring PL, Gaines-Das R, Bangham DR. International Reference Preparation of Human Chorionic Gonadotrophin for Immunoassay: potency estimates in various bioassay and protein binding assay systems. *J Endocrinol* 1980;84:295–310.
- 15) Vaitukaitis J. Human chorionic gonadotropin – a hormone secreted for many reasons. *N Eng J Med* 1979;301:324–6.
- 16) Danzer H, Braunstein GD, et al. Maternal serum human chorionic gonadotropin concentrations and fetal sex predictions. *Fertil Steril* 1980;34:336–40.
- 17) Husa RO. *The Clinical Marker HCG*. New York: Praeger Publishers, 1987.
- 18) Greene MF, de Montserrat MF, Tulchinsky D. In: Tietz NW, editor. *Textbook of Clinical Chemistry*. Philadelphia: Saunders, 1986: 1745.
- 19) Braunstein GD, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976;126:678.
- 20) Babson AL. The IMMULITE Automated Immunoassay System. *J Clin Immunoassay*

1991;14:83-8. 21) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 22) Tietz NW, ed. Clinical Guide to Laboratory Tests. 2nd Ed. Philadelphia: W.B. Saunders, 1990. 23) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 24) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 25) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Intraassay Precision (mIU/mL)

	Mean ¹	SD ²	CV ³	<i>n</i>
1	30.9	1.62	5.2%	20
2	84.6	3.97	4.7%	20
3	492	26.5	5.4%	20
4	3329	120	3.6%	10

Interassay Precision (mIU/mL)

	Mean ¹	SD ²	CV ³	<i>n</i>
1	37	3.2	8.6%	14
2	106	9.5	9.0%	10
3	171	17	9.9%	14
4	3569	279	7.8%	16

Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	52.9	—	—
	4 in 8	26.4	26.5	100%
	2 in 8	13.5	13.2	102%
	1 in 8	6.1	6.6	92%
2	8 in 8	141	—	—
	4 in 8	76.5	70.5	109%
	2 in 8	37.6	35.3	107%
	1 in 8	16.0	17.6	91%
3	8 in 8	598	—	—
	4 in 8	281	299	94%
	2 in 8	155	150	103%
	1 in 8	82.3	74.8	110%
4	8 in 8	3719	—	—
	4 in 8	1976	1860	106%
	2 in 8	936	930	101%
	1 in 8	431	465	93%

Recovery (mIU/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	23.2	—	—
	A	196	186	105%
	B	373	364	102%
	C	798	767	104%
2	—	58	—	—
	A	197	219	90%
	B	379	397	95%
	C	761	800	95%
3	—	318	—	—
	A	475	466	102%
	B	654	644	102%
	C	972	1047	93%

Specificity



Compound ¹	Amount Added ² (ng/mL)	% Cross reactivity ³
FSH	26.8	ND
LH	16.5	ND
TSH	860	ND

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
© 2019-2021 Siemens Healthcare Diagnostics. All rights reserved.
ADJUSTOR AB